




HADER
SOLUTIONS



Implant Buddy

USER GUIDE




WWW.HADER.EU

SERVICE 99% OF IMPLANTS WITH THE HADER IMPLANT BUDDY

Introducing the Implant Buddy, a revolutionary kit that simplifies dental prosthetic work.

This compact and practical solution houses all the essential drivers needed to service 99% of implants from leading brands in the market. Its small size not only fits seamlessly into the dental chair, conserving valuable space, but also optimizes storage in the autoclave. A must-have for dentists and dental technicians, the Implant Buddy streamlines the placement and removal of prosthetic elements over implants, ensuring efficiency and precision in every procedure.



Implant Buddy

Available in three versions


IMPLANT BUDDY DENTAL

Prosthetic connections included:

- SCS Screwdriver – Straumann
- Angled Access Driver
- Unigrip Driver – Nobel Biocare
- Hex 1.3mm – BioHorizons, Astra
- Hex 1.25mm – Zimmer, Camlog
- Hex 1.2mm – 3i Systems, Megagen, BTI
- Hex 0.9mm – Ankylos



Includes a convenient thumb knob for finger pressure



Available references:

1012046-1 - Implant Buddy Dental

1012045 - Implant Buddy Dental with OMEGA torque wrench.

Short Drivers

Total length 17mm
Shaft length 6.5mm

Long Drivers

Total length 24mm
Shaft length 13.5mm

IMPLANT BUDDY LAB

Only for laboratory use.

Prosthetic connections included:

- Straight Driver – Straight Systems
- SCS Screwdriver – Straumann
- Angled Access Driver
- Unigrip Driver – Nobel Biocare
- Hex 1.3mm – BioHorizons, Astra
- Hex 1.2mm – 3i Systems, Megagen, BTI
- Hex 0.9mm – Ankylos



Ergonomic handle ensures a comfortable grip.
Designed for optimal ease of use, in the laboratory



Available references:

1012047-1 - Implant Buddy Dental

Short Drivers

Total length 24mm
Shaft length 16mm

Long Drivers

Total length 40mm
Shaft length 32mm

IMPLANT BUDDY LATCH

Prosthetic connections included:

- Straight Driver – Straight Systems
- SCS Screwdriver – Straumann
- Angled Access Driver
- Unigrip Driver – Nobel Biocare
- Hex 1.3mm – BioHorizons, Astra
- Hex 1.2mm – 3i Systems, Megagen, BTI
- Hex 0.9mm – Ankylos
- Multi-Unit Driver



Includes a latch-to-square adaptor to allow use with square connection torque wrenches



Available references:

1012048-1 - Implant Buddy Latch

1012049 - Implant Buddy Latch with OMEGA torque wrench.

Drivers length 24mm

INDIVIDUAL COMPONENTS



Need a specific component? No problem! The Implant Buddy family offers all individual components for separate purchase. Choose from a range including square and latch drivers, thumb knobs, and adaptors, ensuring you have the exact tools you need for every procedure.

Scan the QR code to download our full catalogue of products

